

Engineered for Enhanced Insertion

INPUT Introducer Sheaths

The Input introducer sheath is engineered for easy insertion and to make your job easier

ENHANCED ENGINEERING INCLUDES

Sit-Up Capability

- > Kink-resistant system provides 60° sit-up comfort for patients with an obturator in place.

Locking Hub

- > Threaded dilator-to-hub interlock prevents dilator backout during insertion, and secure attachment of the sterile sleeve.

Enhanced Tip Geometry

- > Advanced radiofrequency shaping technology creates a long, continuously smooth taper, for minimum resistance during insertions.

Large Internal Geometry

- > Input TS's exceptionally large internal geometries facilitate hemodynamic monitoring and fluid administration.

Precision Engineering

- > Precisely engineered materials and design provide optimal energy transfer while minimizing friction during insertion.

INPUT TS

Engineered to Reduce Drag and Insertion Force

- > FEP sheath with silicone coating enhances the low drag coefficient of the sheath and dilator system.
- > Polyethylene dilator with silicone coating and gradual taper reduces force during insertion.

INPUT PS

Based on the Input TS platform, the Input PS sheath includes a proprietary HydroPel coating that reduces insertion force.

HydroPel Proprietary Hydrophilic coating

- > Hydrophilic coating activates on contact with fluids to produce a smooth, slick, permanent surface.
- > Uncoated near hub to allow for smooth insertion.
- > Pebax material shaft adds the right balance of strength and flexibility for atraumatic insertion into tortuous anatomies.

7 cm sheath

Product Number	French Size	Obturator Included	Sheath Length (cm)	Wire Style
040028A	5	No	7	0.038 DDJ
040031A	6	No	7	0.038 DDJ
040035A	7	No	7	0.038 DDJ

11 cm sheath

Product Number	French Size	Obturator Included	Sheath Length (cm)	Wire Style
051101A	5	No	11	0.038 DDJ
061101A	6	No	11	0.038 DDJ
071101A	7	No	11	0.038 DDJ
081101A	8	No	11	0.038 DDJ
091101A	9	No	11	0.038 DDJ
010101A	10	No	11	0.038 DDJ
011101A	11	No	11	0.038 DDJ

23 cm sheath

Product Number	French Size	Obturator Included	Sheath Length (cm)	Wire Style
060037A	6	Yes	23	N/A
070037A	7	Yes	23	N/A
080037A	8	Yes	23	N/A
090037A	9	Yes	23	N/A

SHEATH WITH NEEDLE AND GUIDEWIRE

11 cm sheath

Product Number	French Size	Sheath Length (cm)	Needle Length (cm)	Needle Gauge
051102A	5	11	7	18
061102A	6	11	7	18
071102A	7	11	7	18
081102A	8	11	7	18
091102A	9	11	7	18

11 cm sheath

Product Number	French Size	Obturator Included	Sheath Length (cm)	Wire Style
050011	5	No	11	0.038 DDJ
060011	6	No	11	0.038 DDJ
070011	7	No	11	0.038 DDJ
080011	8	No	11	0.038 DDJ
090011	9	No	11	0.038 DDJ
853311	5	No	11	N/A
863311	6	No	11	N/A
873311	7	No	11	N/A
883311	8	No	11	N/A
893311	9	No	11	N/A

23 cm sheath

Product Number	French Size	Obturator Included	Sheath Length (cm)	Wire Style
060023	6	Yes	23	N/A
070023	7	Yes	23	N/A
080023	8	Yes	23	N/A
090023	9	Yes	23	N/A

INTRODUCER SHEATH STERILE SLEEVES

Product Number	Sheath Length (cm)
026000	26
060000	60

Medtronic Vascular
3576 Unocal Place
Santa Rosa, CA 95403
USA
Tel: +1.707.525.0111
www.Medtronic.com

Medtronic BV
Earl Bakkenstraat 10
6422 PJ Heerlen
The Netherlands
Tel: +31.45.566.8000
Fax: +31.45.566.8668

Product Services
Tel: +1.888.283.7868
Fax: +1.800.838.3103

*Check with sales representative for availability. Prior to use, please consult labels and inserts for any indications, contraindications, hazards, warnings, cautions, and use.

CAUTION: Federal Law restricts this device to sale by or on the order of a licensed healthcare practitioner.

© 2005 Medtronic, Inc. All rights reserved. Printed in USA. UC200600929EN 10/05

INPUT®
Introducer Sheaths



Engineered for Enhanced Insertion

ACCESSORIES

OBTURATOR

- > Helps maintain sheath patency and protect from outside contamination.
- > Twist'N Shut locking system locks obturator in place to prevent backout from sheath.

VESSEL DILATOR

- > Polyethylene dilator for use in arterial or venous applications.

STERILE SLEEVE

- > Maintains sterility of indwelling catheter shaft proximal to the introducer hub.
- > Twist'N Shut locking system locks sleeve into input hub.