




Coronary Artery Bypass Surgery "en miniature"


ValveGate™ Mini



Efficient surgery. Less invasive. Optimum results.



ValveGate™ Mini



A

B

C

D

E

F

G

H

I



J K L M N O

- | | | |
|---|----------|--|
| A | 34-7820S | ValveGate™ DeBakey Grasper , Ø 5 mm
1x2 narrow 1.5x11 mm, 17 cm |
| B | 34-7821S | ValveGate™ DeBakey Grasper angled upwards
Ø 5 mm, 1x2 narrow 1.5x14mm, 17cm |
| C | 34-7822S | ValveGate™ DeBakey Grasper angled sideways
Ø 5 mm, 1x2 narrow 1.5x14 mm, 17 cm |
| D | 34-7825S | ValveGate™ DeBakey Grasper , Ø 5 mm
1x2 wide 2.8x11 mm, 17 cm |
| E | 34-7827S | ValveGate™ DeBakey Grasper , Ø 5/3 mm
malleable shaft, 1x2 narrow 1.5x11 mm, 17 cm |
| F | 34-7826S | ValveGate™ DeBakey Grasper soft, Ø 5 mm
2x3 extra wide, 4.8x14 mm, 17 cm |
| G | 34-7875S | ValveGate™ Mini DeBakey Grasper , Ø 3.4 mm
1x2 1.5x15 mm angled up, 17 cm |
| H | 34-7876S | ValveGate™ Mini DeBakey Grasper , Ø 3.4 mm
1x2 1.5x15 mm straight, 17 cm |
| I | 34-7881S | ValveGate™ DeBakey Grasper , Ø 5 mm
1x2 narrow 1.5x15 mm straight, 17 cm |
| J | 34-7800S | ValveGate™ TC Needle Holder straight, Ø 5 mm
with ratchet, 17 cm |
| K | 34-7802S | ValveGate™ TC Needle Holder curved, Ø 5 mm
with ratchet, 17 cm |
| L | 34-7805S | ValveGate™ TC Ryder Needle Holder straight, Ø 5 mm
1.2 mm tip, with ratchet, 17 cm |
| M | 34-7806S | ValveGate™ Mini TC Needle Holder curved, Ø 2.7 mm
with ratchet, 17 cm |
| N | 34-7804S | ValveGate™ TC Ryder Needle Holder straight, Ø 3 mm
0.8 mm tip, with ratchet, 17 cm |
| O | 34-7808S | ValveGate™ Mini TC Needle Holder , angled upwards
Ø 2.7 mm, with ratchet, 17 cm |




ValveGate™ Mini







- | | | |
|---|----------|--|
| A | 34-7830S | ValveGate™ Scissors curved 15°, Ø 5 mm, 17 cm |
| B | 34-7835S | ValveGate™ Scissors curved 30°, Ø 5 mm, 17 cm |
| C | 34-7840S | ValveGate™ Scissors curved 70°, Ø 5 mm, 17 cm |
| D | 34-7860S | ValveGate™ Potts Scissors 125°, Ø 3.4 mm
6 mm blade, 17 cm |
| E | 34-7862S | ValveGate™ Potts Scissors 45°, Ø 3.4 mm
6 mm blade, 17 cm |
| F | 34-7863S | ValveGate™ Mini Potts Scissors 125°, Ø 3 mm
4 mm blade, 17 cm |
| G | 34-7864S | ValveGate™ Mini Potts Scissors 45°, Ø 3 mm
4 mm blade, 17 cm |
| H | 34-7865S | ValveGate™ Mini Dissecting Scissors slightly curved
Ø 3 mm, 17 cm |
| I | 34-7850S | ValveGate™ Knot Pusher, Ø 5 mm
17 cm, push to close |
| | 34-7849S | ValveGate™ Knot Pusher, Ø 5 mm
17 cm, push to close, short jaw |
| J | 34-7851S | ValveGate™ Knot Pusher, Ø 5 mm
17 cm, push to open, short jaw |
| K | 34-7852S | ValveGate™ Mini Knot Pusher, Ø 3mm
push to close, 17 cm |
| L | 03-7611 | Chest Tube Passer 36 cm with ratchet
(used as retractor to depress lung), with ring handle |
| M | 34-7878S | ValveGate™ Mini IMA Holding Forceps upwards angled
0.8 mm jaw with 1x2 mouse teeth, Sapphire™ jaws,
17 cm, Ø 3.4 mm |
| N | 34-7829S | ValveGate™ Ring Tissue Forceps, Ø 5 mm
Sapphire™ jaws, 17 cm |
| O | 34-7460S | ValveGate™ Clip Applying Forceps, Ø 5 mm
for Vitalitec/SLS clip small, 17 cm |

Alternative lengths B=12 cm and C=15 cm instead of S=17 cm available at same price.
12/15 cm might be beneficial for CABG.



We are working in the three business fields „Cardio“, „Neuro“ and „Powered“. In all these fields we want to offer our valued clients a complete range of products and services. Our products are clearly described and illustrated in well-structured sales documents. A small selection of various catalogs we show you below.







GEISTER Medizintechnik GmbH

Föhrenstrasse 2 | D-78532 Tuttlingen | Germany

☎ +49 7461/966240 | 📠 +49 7461/9662422

www.geister.com | info@geister.com